# *Amnio* slide

**Protocollo di utilizzo**  
**Protocole d'utilisation**  
**Benutzungsprotokoll**  
**Protocol of use**  
**Protocolo de utilización**  
**Protocolo de utilização**  
**Πρωτόκολλο χρήσης**

- PER COLTURADI CELLULE
- POUR LA CULTURE DES CELLULES
- FÜR ZELLKULTUREN
- FOR CELL CULTURE
- PARA CULTIVO DE CÉLULAS
- PARA CULTURA DE CÉLULAS
- ΓΙΑ ΚΑΛΛΙΕΡΓΕΙΑ ΚΥΤΤΑΡΩΝ

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## Symbols used in the labelling

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English

## Chamber for culture and analysis of amniocytes and chorionic villi

Ref. EK AMS 60F SuperFrost®, 60 tests

### Material supplied

Each case contains 30 ready to use trays, 60 chambers, treated with Ethylene Oxide.  
Validity: 2 years.

### Amniocyte culture in situ protocol

Note: thaw the Amniomed Plus (ref. EKAMG200) in a 37°C water bath and mix well by swirling prior to use. Addition of L-Glutamine, antibiotics and serum are not necessary since these components are already present.

1. Place 1 ml of cell suspension gently onto each Amnioslide, on the microscope slide.
2. Incubate a 37°C, 5% CO<sub>2</sub> for 6 - 12 hours, with minimal mechanical disturbance.
3. Add 4 ml of fresh Amniomed Plus.
4. On Day 6, replace the culture medium with 5 ml of fresh Amniomed Plus.
5. On Day 7 / Day 8 check for progress of growth using an inverted microscope. When the areas of growth show a number of cells in mitosis equal to three to four microscope fields (100X), add 30 µl of Colcemid (10 µg/ml).
6. Incubate at 37°C, 5% CO<sub>2</sub> for 60 - 90 minutes.

### Chromosome preparation, manual method

1. Suck the medium up completely.
2. Drop 6 ml of hypotonic solution (distilled H<sub>2</sub>O, 0,6% sodium citrate, 0,1% KCl) and incubate at room temperature for 10 minutes.
3. Suck the hypotonic solution up and drop 6 ml of Ibraimov solution (distilled H<sub>2</sub>O, 5% acetic acid), incubate at room temperature for 5 minutes.
4. Suck the Ibraimov modified solution up and drop 6 ml of fresh fixative mixture (methanol : acetic acid, 4 : 1) at room temperature. The time for this wash has no influence on final results.
5. Repeat 4) twice.
6. Remove the microscope slide and dry under constant condition of right temperature and relative humidity into Optichrome (ref. EKAMH950).

### Chromosome preparation, automated method

1. Put Amnioslide in Autochrome (ref. EKEKAMH1001).
2. Chose the proper program.
3. At the end of the automated process remove the microscope slide and dry under constant condition of temperature and relative humidity into Optichrome.

### Chorionic villi sampling (CVS) culture protocol

1. Transfer the specimen into a 60 mm Petri dish containing 5 ml of RPMI 1640.
2. Wash villi with fresh medium to remove blood cells.
3. Using the inverted microscope, carefully dissect any remaining clots or decidual fragments.

### In situ culture

1. Transfer the villi in a 15 ml sterile centrifuge tube containing 1 ml of Pronase E (4\*10<sup>6</sup> PU/g) and incubate at room temperature for 4 - 6 minutes.
2. Add 3 - 5 ml of cold Hank's Balanced Salts Solution (kept at +4°C).
3. Centrifuge for 5 minutes at 1500 rpm and discard the supernatant.
4. Add 2 ml of sterile Collagenase type II (1 mg/ml) and incubate at 37°C for 10 minutes.
5. Add 3 - 5 ml of cold HBSS (kept at +4°C).
6. Centrifuge for 5 minutes at 1500 rpm and discard the supernatant.
7. Add 2 ml of Amniomed Plus and resuspend the pellet.
8. Prepare 2 to 6 Amnioslide, depending on the number of cells. Add 4,5 ml of Amniomed Plus medium to each Amnioslide.
9. Drop about 500 µl of the cell suspension on each Amnioslide.
10. Incubate at 37°C in 5% CO<sub>2</sub>.
11. On Day 5 replace the culture medium with 5 ml of fresh Amniomed Plus.
12. On Day 6 / Day 7 check for progress of growth using an inverted microscope. When the areas of growth show a number of cells in mitosis equal to three to four microscope fields (100X), add 30 µl of Colcemid (10 µg/ml).
13. Incubate at 37 °C, 5% CO<sub>2</sub> for 4 - 6 hours.

### Chromosome preparation, manual or automated method

See protocol chromosome preparation of Amniotic fluid.

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**EuroClone S.p.A.**