Cell Culture Slide

1 well and 2 well

REF 101307, 101308

INTENDED USE

In situ culture of amniocytes and chorionic villi derived cells for prenatal cytogenetic diagnosis.

INTRODUCTION

- SPL Cell Culture Slide is a polystyrene chamber with cover, assembled to a borosilicate microscope glass slide by a slide holder.
- At the end of insitu culture, the chamber is removed from the slide for further necessary experiment.
- To determine chromosomal or genetic disorders in the fetus, some prenatal cytogenetic techniques such as chorionic villus sampling or aminocentesis have been performed.

MATERIAL SUPPLIED

- Each box contains 2 sterile packages ready for use, for a total of 12 cell culture slide.
- Slide dimension: 76 x 26 mm
- Culture area : 11.2cm² for 1 well, 10.4cm² for 2 well
- Validity: 3 years

MATERIAL SUPPLIED BUT NOT PROVIDED

- Cell culture media for amniocytes and chorionic villi derived cells
- Adjustable, automatic micropipettes with disposable tips

PROTOCOL CELL CULTURE SLIDE

English

Code 101307 - CULTURE SLIDE 1 WELL

Amniocytes- Culture preparation

1. Centrifuge the amniotic fluid sample at 2000 rpm for 10 min, discard the supernatant, leaving 2 ml of liquid and resuspend the pellet. Place 0.5 ml of cell suspension in the middle of the well, and add 2-2.5 ml of fresh culture medium.
2. Incubate at 37°C , 5% CO²
3. At day 6° completely change the culture medium with 2.5-3 ml of fresh culture medium, and incubate.
4. On day 7° check the growth through an inverted microscope. When the growth areas show a sufficient number of cells in mitosis add Colcemid to a final concentration of 0.05 ug/ml .
5. Incubate at 37°C , 5% CO² for about 3 hours.

Chromosome fixation and slide preparation

1. Remove completely the medium and remove the slide from the chamber.
2. Dip the slides in a Coplin jar containing: 40 ml of hypotonic solution 1% Sodium Citrate and 40 ml of hypotonic solution 0.3% NaCl; incubate at room temperature for 10 minutes.
3. Dip the slides in a Coplin jar containing 80 ml of Ibraimov solution (distilled H₂O, 5% acetic acid), incubate at room temperature for 10 minutes.
4. Dip the slides in a Coplin jar containing 80 ml of fixative solution (methanol:acetic acid, 3:1) at room temperature.
The time of this incubation does not influence the final result.
5. Repeat step 4 twice. For the last wash the fixative must be prepared shortly before use.
6. Proceed and adequate drying in constant conditions of temperature and relative humidity (25°C , 35% RH).

CHORIONIC VILLI (CVS) – Culture preparation

1. Transfer the specimen into a 60 mm Petri dish containing 5 ml of Hank's (HBSS) without Ca²⁺-Mg²⁺.
2. Wash villi with fresh medium to remove blood cells.
3. Using the inverted microscope, carefully dissect any remaining clots or decidual fragments.
4. Transfer the villi in a 15 ml sterile centrifuge tube and centrifuge for 2 minutes at 2000 rpm; discard the supernatant and add 1 ml of Pronase E(4*106PU/g); incubate at room temperature for 10 minutes on a shaker.
5. Add 5 ml of cold Hank's (HBSS) without Ca²⁺-Mg²⁺(kept at +4°C).
6. Centrifuge for 5 minutes at 2000 rpm and discard the supernatant.
7. Add 2 ml of Collagenase type II (1 mg/ml) and incubate at 37°C for 10 minutes.
8. Add 3-5 ml of cold (HBSS) without Ca²⁺-Mg²⁺(kept at +4°C).
9. Centrifuge for 5 minutes at 2000 rpm and discard the supernatant.
10. Repeat points 8 and 9.
11. Add 1 ml of complete medium for villi and resuspend the pellet.

In situ culture

1. Prepare the Cell Culture Slide, seed 5 drop of the cell suspension on each Slide and check the concentration of the cells. Add 2,5-3 ml of complete medium for villi to each Slide.
2. Incubate at 37°C in 5% CO²
3. On Day 6 replace the culture medium with 2,5-3 ml of complete medium for villi.
4. On Day 7 check for progress of growth using an inverted microscope.
At the optimal cell growth, add Colcemid (10 µg/ml) at the final concentration of 0,05ug/ml.
5. Incubate at 37 °C, 5% CO² for 4 - 6 hours.

Follow the protocol of "Chromosome fixation and slide preparation" for Amniocytes.

Code101308- CULTURE SLIDE 2 WELLS

Amniocytes- Culture preparation

1. Centrifuge the amniotic fluid sample at 2000 rpm for 10 min, discard the supernatant, leaving 2 ml of liquid and resuspend the pellet. Place 0.5 ml of cell suspension in the middle of the well, and add 1,5 ml of fresh culture medium.
2. Incubate at 37 ° C , 5 % CO²
3. At day 6° completely change the culture medium with 1,5 ml of fresh culture medium, and incubate.
4. On day 7° check the growth through an inverted microscope. When the growth areas show a sufficient number of cells in mitosis add Colcemid to a final concentration of 0.05 ug/ml .
5. Incubate at 37 ° C , 5 % CO² for about 3 hours.

Chromosome fixation and slide preparation

1. Remove completely the medium and remove the slide from the chamber.
2. Dip the slides in a Coplin jar containing: 40 ml of hypotonic solution 1% Sodium Citrate and 40 ml of hypotonic solution 0.3% NaCl; incubate at room temperature for 10 minutes.
3. Dip the slides in a Coplin jar containing 80 ml of Ibraimov solution (distilled H₂O, 5% acetic acid), incubate at room temperature for 10 minutes.
4. Dip the slides in a Coplin jar containing 80 ml of fixative solution (methanol:acetic acid, 3:1) at room temperature.
The time of this incubation does not influence the final result.
5. Repeat step 4 twice. For the last wash the fixative must be prepared shortly before use.
6. Proceed and adequate drying in constant conditions of temperature and relative humidity (25°C , 35% RH).

CHORIONIC VILLI (CVS) – Culture preparation

1. Transfer the specimen into a 60 mm Petri dish containing 5 ml of Hank's (HBSS) without Ca²⁺-Mg²⁺.
2. Wash villi with fresh medium to remove blood cells.
3. Using the inverted microscope, carefully dissect any remaining clots or decidual fragments.
4. Transfer the villi in a 15 ml sterile centrifuge tube and centrifuge for 2 minutes at 2000 rpm; discard the supernatant and add 1ml of Pronase E (4*106PU/g); incubate at room temperature for 10 minutes on a shaker.
5. Add 5 ml of cold Hank's (HBSS) without Ca²⁺-Mg²⁺ (kept at +4°C).
6. Centrifuge for 5 minutes at 2000 rpm and discard the supernatant.
7. Add 2 ml of Collagenase type II (1 mg/ml) and incubate at 37°C for 10 minutes.
8. Add 3-5 ml of cold (HBSS) without Ca²⁺-Mg²⁺ (kept at +4°C).
9. Centrifuge for 5 minutes at 2000 rpm and discard the supernatant.
10. Repeat points 8 and 9.
11. Add 1 ml of complete medium for villi and resuspend the pellet.

In situ culture

1. Prepare the Cell Culture Slide, seed 3 drops of the cell suspension on each Slide and check the concentration of the cells. Add 1,5 ml of complete medium for villi to each Slide.
2. Incubate at 37°C in 5% CO²
3. On Day 6 replace the culture medium with 1,5 ml of complete medium for villi.
4. On Day 7 check for progress of growth using an inverted microscope. At the optimal cell growth, add Colcemid (10 µg/ml) at the final concentration of 0,05ug/ml.
5. Incubate at 37 °C, 5% CO² for 4 - 6 hours.

Follow the protocol of "Chromosome fixation and slide preparation" for Amniocytes.

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SPL Life Sciences Co., Ltd 26, Geumgang-ro 2047 beon-gil, Naecheon-Myeon, Pocheon-si, Gyeonggi-do 487 835, Korea
EUROCLONE S.p.A. Via Figino, 20/22 - 20016 Pero (MI) Italy