Chromosome *kit* Φ

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Protocol of use

Protocolo de utilización

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Complete medium for culture of peripheral blood cells

Material supplied

EKAMTP, 10 x 5 ml culture tubes
EK AMTB 100, 1 bottle 100 ml
EK AMTB 500, 1 bottle 500 ml

Product stability

Chromosome Kit and Medium P is shipped at -18/-22°C, and has to be stored at -18/-22°C until the expiration date indicated on the product label. Before use, thaw the Medium in a refrigerator at +2°C/+8°C and mix well by swirling prior to use. Once thawed, it can be stored at +2°C/+8°C for maximum 14 days. Chromosome P can be thawed and aseptically transferred into smaller aliquots for convenience. Aliquots can be frozen and thawed at time of use; do not repeat "thaw-freeze" cycles more than twice.

English

Protocol of use

Culture preparation

1. Aseptically add up to 0.5 ml of whole blood to one test tube or dispense 5 ml of bottled medium in a T-Flask and add the sample.
2. Replace the cap of the culture tube.
3. Mix the culture tube (or the T-Flask) contents and incubate at 37°C placing the tube horizontally (on its flat surface). Addition of CO₂ is not necessary as the medium composition is designed for optimal lymphocyte proliferation in a closed system.

Chromosome fixation and slides preparation

1. After 72 hours of culture, add 50 µl of Colcemid (10 µg/ml) and incubate at 37°C for 90 minutes (approx.). This last incubation time may be reduced to 15 minutes if the aim is obtaining a higher number of prometaphases.
2. After incubation, transfer the culture in a tube, if necessary, and centrifuge for 3 - 4 minutes at 2.000 rpm (higher speeds will not damage the cells).
3. Discard the supernatant, taking care not to lose part of the cell pellet (it is recommended to leave approximately 0.5 ml of material in the tube).
4. Vigorously resuspend the pellet.
5. Add 5 ml of hypotonic solution (distilled H₂O, 0,075M KCl) and mix adequately.
6. Incubate for at least 7 minutes at room temperature.
7. Centrifuge at 2.000 rpm for 3 - 4 minutes.
8. Discard the supernatant according to step 3.
9. Vigorously resuspend the pellet.
10. Add 5 ml of Ibraimov solution (distilled H₂O, 5% acetic acid) and mix adequately.
11. Immediately centrifuge at 2.000 rpm for 3 - 4 minutes.
12. Discard the supernatant according to step 3.
13. Vigorously resuspend the pellet.
14. Add 5 ml of fixative solution (methanol or ethanol: acetic acid, 3 : 1) and mix adequately.
15. Repeat steps 11-14.
16. Immediately centrifuge at 2.000 rpm for 3 - 4 minutes.
17. Carefully discard almost all the supernatant.
18. Add a few drops of freshly prepared fixative, taking into consideration that the number of drops should be appropriate to the cell concentration in the pellet.
19. Drop a small aliquot of the cell suspension onto a clean and wet microscope slide.
20. Evaporate under constant condition of right temperature and relative humidity into Optichrome PLUS (cod. EKAMH960).

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EuroClone S.p.A.