


ChromoQuant[®] AZF

PCR/CAPILLARY ELECTROPHORESIS FOR CHROMOSOME Y MICRODELETIONS MOLECULAR DIAGNOSTIC

AZFa
AZFb (P5/proximal P1)
AZFbc (P5/distal P1 or P4/distal P1)
AZFc (b2/b4)

- ❖ Total of **17 markers** (terminal deletion included)
- ❖ 8 extra markers included for proximal/distal AZFa/AZFb breakpoints
- ❖ Validated Instruments : ABI 310, 3100, 3130, 3500, 3730


CE marked
Directive 98/79/EC
ISO 13485 / ISO 9001



In accordance with EAA
and EMQN Guidelines 2013

Guidelines : European Academy of
Andrology (EAA) / European Molecular
Genetics Quality Network (EMQN)

ChromoQuant[®] AZFa proximal and distal deletion



ChromoQuant AZF Kit (17 markers)

Reference 800.001-32

1 tube 32 tests
7 non polymorphic markers AZFa, b, c
+ 2 controls

1 tube 10 tests
4 Extra-markers AZFa + 2 controls

1 tube 10 tests
4 Extra-markers AZFb + 2 controls